

Origin ships with a collection of examples to show LabVIEW users how to create VIs to communicate between LabVIEW and Origin

The examples provided in the following Origin installation folder:

<exe>\Samples\Automation Server\LabVIEW\

Notes: These examples require Origin 8 SR4 or later builds.

Examples included in this section:

- ◇ Origin Column Set Get Data
- ◇ Batch Linear Regression
- ◇ Simple NLFit
- ◇ On The Fly Analysis
- ◇ Simulate DAQ
- ◇ 6009 Simple LV71
- ◇ Realtime Curve Fitting
- ◇ 6009 LV71
- ◇ Plotting Data
- ◇ Sending XY Data To Origin
- ◇ Simulated Peak
- ◇ Matrix Access
- ◇ Sending Waveforms to Origin
- ◇ Analysis Template
- ◇ Waveform Envelope